

## Stands

- Solid cast base, coated with chemical resistant epoxy
- Largest base accommodates hotplate or large vessels; Leveling knobs provide added adjustability
- Compact base maximizes bench space and fits into smaller enclosed work stations or hoods


## Clamps

- Securely hold a large portion of the stirrer support rod maximizing stability
- Multipurpose - also used with chain or extension clamps to hold accessories



2

A110	A112, A113, A114	A210	Customize
------	------------------	------	-----------



Use with all stirrers except 250  
Accommodates maximum vessel diameters of 323 (12 3/4)

Use with all stirrers except 250  
Choice of three stand rod lengths  
Includes two extra threaded holes to accommodate A300 - 710 x 16 (28 x 5/8) rod for accessories

Use with 2002 and 2010  
Unique offset design centers impeller in vessel  
Accommodates maximum vessel diameters of 190 (7 1/2)

A110BASE - base only  
ROD28  
710 x 25 (28 x 1)  
ROD38  
960 x 25 (38 x 1)  
ROD48  
1200 x 25 (48 x 1)



3



Clamps	A120	A122
	Use with 12.5 (1/2) support rod	Use with 16 (5/8) support rod
Adjustable to accept 15 - 30 ( 9/16 - 1 3/16 ) diameter stand rods		

mm (in)