

# ELECTROLAB

## Bottle Rotating Apparatus



**BETATEK** INC.  
*Instrumentation Superbly Supported*  
*L'instrumentation à son meilleur*

654 Petrolia Road,  
Toronto, ON M3J 2W3  
Tel: (416) 736-6166  
Toll-Free: (800) 387 3570  
Fax: (416) 736-9346


## Dissolution Testers



### Bottle Rotating Apparatus (ERB-16W)

#### Water Bath Heating System

CE



**ERB-40W**  
(supporting 40 bottles of 10 ml each)  
is also available

#### Key Features

- Supports up to 16 bottles of 100 ml each
- Magnetically coupled drive and precise speed control from 5 to 50 RPM with soft start facility
- Programmable temperature in the water bath from 30°C to 50°C
- Snap fitting of bottles for easy placement and removal
- Robustly built clear polycarbonate water bath to withstand high temperatures

### Air bath Heating System (ERB-40A)



**ERB-16A**  
(supporting 16 bottles of 100 ml each)  
is also available

#### Key Features

- Supports 40 bottles of 10 ml each
- Precise speed control from 5 to 50 RPM with soft start facility
- Forced air heating system with programmable temperature range from 30°C to 50°C
- 24 x 7 continuous operation
- Snap fitting of bottles for easy placement and removal
- Impact resistant clear polycarbonate air bath to withstand high temperatures