(MODEL_{2500rtd}

The Distek Model 2500 RTD Dissolution System offers great flexibility and configurability. It can be configured as USP Apparatus 1, 2, 5 and 6 plus intrinsic dissolution. The extensive users, method and reporting storage capabilities are enabled by a brilliant touch screen. Patented wireless temperature sensors continuously monitor and display the in-vessel temperature for each vessel. Unparalleled control of the instrument plus automatic storage of component serial numbers and automatic qualification reminders make the Model 2500 RTD robust and reliable enough for QC applications yet flexible enough for R&D use.

WIRELESSIN-SHAFT TEMPERATURE SENSORS

Patented wireless temperature sensors control and monitor the temperature.

- Continuously monitor, display and record the in-vessel temperature for each vessel
- Ensure and document temperature compliance throughout the entire dissolution test
- Eliminate additional labor to measure actual temperature in vessels

VOLUME FLEXIBILITY

Accommodates up to 8 vesselsranging in size from 100 mL (utilizing Distek's small volume conversion kit) to 2L, allowing the user to run a diverse range of tests.

• SELF-PRIMING CIRCULATOR/ HEATER

The self-priming design and precise microprocessor control allow for improved temperature uniformity and reduced heating times. The low-profile design allows the circulator to be stowed behind the bath, requiring no additional bench space.

ELECTRONICQUALIFICATION GUIDANCE

Settable automatic electronic instrument requalification reminders eliminate missed validation dates and ensures continuous instrument uptime.

COLOR TOUCH SCREEN DISPLAY

The icon driven user interface lowers overall cost by reducing training time and user errors while maximizing productivity and command of the dissolution test.

• QUICK CHANGE PADDLESAND BASKETS

Interchangeablestirring elements coupled with the automatic height adjustment provide a fast transition between Apparatus 1 and 2 without having to remove the shaft.

AUTOMATIC HEIGHT ADJUSTMENT

Setting the height is done through a one-time adjustment and utilization of the interchangeable paddles and baskets

• WAKE UP AND SLEEPMODE

Programmable wake up and sleep modes allow the user to schedule thermocirculator start up and shutdown times to save time and conserve energy.







SPECIFICATIONS

Dissolution Vessels	Six VesselsStandard (7 or 8 Optional)
Volume	25 - 2000 mL (Volumes below 500 mL Require Low Volume Kit)
Water Bath	One Piece Molded Acrylic
RPM Control Range	25 - 300 RPM, Digitally Controlled, Closed Loop
RPM Resolution	0.1 RPM
RPM Accuracy	± 1 RPM up to 100 RPM and $\pm 1\% > 100$ RPM
Motor	High Torque, Permanent Magnet
Vessel Temperature Control	Ambient to 55°C, TCS-0500Thermocirculator
Temperature Resolution	0.01°C
Temperature Accuracy	±0.25°C
Shaft Wobble	Lessthan 0.010" / 0.254mm (Total Indicator Runout)
Lab Temperature Min/Max	20°C / 25°C
Lab Humidity Min/Max	20%/ 80% Relative Humidity
Display	5.7" Color TouchScreen
User Management	Manage up to 50 Userswith Multiple AccessLevels
	• Manual
Program Modes	Automatic (Up to 100 Pre-Programmed Methods)
	Distek EVO 4300, Eclipse 5300, Opt-Diss 405, Opt-Diss 410
Interface Ports	USB, Ethernet, RS-232, RS-485
Construction Materials	Cast Aluminum, Stainless Steel, Engineered Plastic
Dimensions	27.5" (w) x 38.5" (h) x 19.5" (d) / (70 cm x 98 cm x 50 cm) TCS: 11.5" (w) x 8.5" (h) x 7" (d) / (29 cm x 22 cm x 18 cm)
Weight	80 lbs. / (36 kg) TCS: 7 lbs. / (3 kg)
Electrical Power	115V ± 15V 50/60 Hz 10A or 230V ± 15V 50/60 Hz 6A (Operating Voltage Pre-Set at Factory) TCS: 115V ± 15V 50/60 Hz 10A or 230V ± 15V 50/60 Hz 7A



Dashboard screen when wireless in-shaft temperature sensorsare Disabled



General Settings screen for both Enabled & Disabled in-shaft temperature sensors



Dashboard screen when wireless in-shaft temperature sensors are Enabled





654 Petrolia Road, Toronto, ON M3J 2W3 Tel: (416) 736-6166 Toll-Free: (800) 387 3570 Fax: (416) 736-9346

