CIPHER SOFTWARE

The Cipher dissolution control software is a Windows based application designed and developed by Distek which can be installed locally or networked providing full access to any Distek dissolution instrument on the network. Cipher offers the flexibility to set the level of security appropriate for your laboratory. When enabling all available security features and audit trails and combined with the appropriate customer developed policies, Cipher will meet 21 CFR part 11 compliance. Easy to use and intuitive, the Cipher software features make preparing, running, monitoring, and reporting dissolution test results quick and easy.

21 CFR PART 11 COMPLIANCE

Combined with the appropriate customer policies, Cipher balances the requirements of compliance with limited intrusion on user activities.

METHOD WIZARD

The user is walked through method setup for a dissolution instrument and autosampler.

• LIMS COMPATIBLE

When configuring the system, two clicks of the mouse permit the automatic exporting of result files in an ASCII format for access and use by a LIMS package.

• eSIGNATURES

The use of electronic signatures guarantees that only approved methods can be used and reports released.

• INSTRUMENT NETWORKING

Users are able to log on to a client work station and view the current status of any Distek system on the network. Permits real time, remote instrument setup and monitoring.

AUDIT TRAIL

Event, system and method changes are tracked with a computer generated date and time stamp. Additional information logged includes the user, old and new values, and the reason for change.

• COMPENDIA ALERTS

The user can set the specific tolerance parameters for rpm and temperature requiring the system to alert the user should a reading fall outside the acceptable range.







CIPHER

SPECIFICATIONS

Window Requirements

For Windows 7 / 10, you must have administrative privileges to install Distek Cipher Software.

Full Administrative access must also be available to the network server (domain controller) if the application will be installed in a Server/ Client configuration. Access to the Server User Manager is required for adding Groups and granting User access to those Groups.

Minimum Computer Requirements (Client): • Intel Pentium ® Processor 1.8 GHz or equivalent • Microsoft ® Windows 7 / 10 (32 or 64 bit) • Microsoft .NET Framework 4.5.1 or higher • Internet Explorer 11 or higher • Minimum of 2GB of RAM • 200MB of available hard disk space • Available USB port (for a serial hub)		Minimum Server Requirements (Only required for a Server): • Intel Pentium ® Processor 1.8 GHz or equivalent • Microsoft ® Windows Compatible Server • Microsoft .NET Framework 4.5.1 or higher • Internet Explorer 11 or higher • Minimum of 4GB of RAM • 200MB of available hard disk space	
TCS Model 2100 C Model 2500 / 2500 RTD Model 2500 Select Premier 5100 Evolution 6100 Evolution 6300 symphony 7100	FW Rev. 1.40 or higher FW Rev. 2.31 or higher FW Rev. 2.05 or higher FW Rev. 3.00 or higher FW Rev. 3.00 or higher FW Rev. 2.12 or higher FW Rev. 2.12 or higher FW Rev. 2.02 or higher	DS 4300 Eclipse 5300 Evolution 4300	FW Rev. 2.10 or higher FW Rev. 1.02 or higher FW Rev. 2.10 or higher
Distek Cli	Distek	VORK SERVER Security Server	Distek Client &
Instrument (LAB #	t Server		Instrument Server (LAB #2)



654 Petrolia Road, Toronto, ON M3J 2W3 Tel: (416) 736-6166 Toll-Free: (800) 387 3570 Fax: (416) 736-9346

SGS