

Enzyme Indicator Technology

An Exclusive Partnership with



TOMI[™] now offers Protak's line of Enzymatic Biological Indicator (EI) products. This validation option is now available for SteraMist[™] decontamination projects both in the U.S. and abroad, with **TOMI[™] being the only supplier in the U.S.** TOMI[™] has the opportunity to offer its clients not simply the most innovative technology on the market for decontaminating spaces, but now the fastest validation process to prove it.

About Protak Enzyme Indicator Technology:

- Protak is the globally exclusive appointed licensee of thermostable Adenylate Kinase technology, developed over 15 years by Public Health England, an executive agency of the United Kingdom's Department of Health.
- The only manufacturer of El's for hydrogen peroxide decontamination performance validation in the world.
- El technology offers the pharmaceutical and life sciences market a **60 second quantitative validation process post decontamination treatment** vs. the standard in today's market which is a binary pass/fail with a 7 day wait prior to reusing the treated room or releasing the sterilized products.
- The El package includes: El's (12pk or 192 pk), ADP Reagent Set, choice of Twin Shot Carousel Reader or 160 Chain Reader, and a Flight Case.



info@tomimist.com 800.525.1698 www.tomimist.com





Benefits of Enzyme Indicators (El's):

- Fast Results in Less than 60 seconds: Instant results + instant action = what's needed to reduce waste, improve productivity, and save lives.
- Comprehensive & Accurate Data Capture: El readers provide quantifiable, linear readings that provide users with substantially more information than the current Biological Indicators (BI's).
 El's can act as early warning systems

 allowing performance monitoring, threshold release and automatic data capture.
 - **Performance Assurance:** El's are manufactured to eliminate rogue results and provide more accurate data than conventional BI tests. El test packages also include an optional on-site self-validation process.
 - **Simple & Easy to Use:** A simple test strip is used alongside an luminometer, which reads and delivers results instantly.
 - **Highly Cost Effective:** The many benefits of the initial investment in El reader and strips include immediate results, longer and stable shelf life of the indicators, and the ability to resume same day operations.
 - **Proven & Researched:** Created and Extensively tested by Public Health England (PHE) for 15 years, El Technology is protected globally with five families of over 60 patents.

654 Petrolia Road, Toronto, ON M3J 2W3

Tel: (416) 736-6166 Toll-Free: (800) 387 3570